

# INTRODUCTION TO FLUORESCENCE MICROSCOPY


- Emmanuel MARGEAT, Research Director CNRS
- Centre de Biochimie Structurale
- Office B4, Building 29
- [margeat@cbs.cnrs.fr](mailto:margeat@cbs.cnrs.fr)
- *Teacher for Intro to qbio, Imaging Biological systems*
- *Responsible for Lab1*


- Draw and describe a Jablonsky-Perrin diagram
- Why is the fluorescence emission wavelength always larger than the excitation wavelength? What is the Stokes shift?
- Notions about the excited state lifetime and how to measure it
- Pro and cons of fluorescent proteins and dyes.
- What is the maximum resolution attainable in a widefield microscope? Calculate in terms of NA and wavelength.
- Draw the optical path of a fluorescence microscope. Describe the components

## JABLONSKI – PERRIN DIAGRAM

- Singlet States & Triplet State
- Vibrational states
- Vibrational relaxation
- Internal conversion
- Intersystem crossing
- Radiative and non-radiative decays & rates
- Absorption & timescale
- Fluorescence & timescale
- Phosphorescence & timescale
- Quantum yield equation
- Excited state lifetime equation
- Origin of the excitation spectrum
- Origin of the emission spectrum
- The Stokes shift
- The case of FRET
  - FRET transition & rate
  - Lifetime equation in the presence of FRET
  - FRET equations

## JABLONSKI – PERRIN DIAGRAM







QUANTUM  
YIELD

$$Q = \frac{k_r}{k_r + k_{nr}}$$

EXCITED STATE  
LIFETIME

$$\tau = \frac{1}{k_r + k_{nr}}$$

## Jablonski diagram in the presence of FRET




$$\tau_D = \frac{1}{k_r + k_{nr}}$$

$$\tau_{D,A} = \frac{1}{k_r + k_{nr} + k_{FRET}}$$


$$E = \frac{k_{FRET}}{k_r + k_{nr} + k_{FRET}} = 1 - \frac{\tau_{D,A}}{\tau_D}$$

# Two techniques for excited state lifetimes measurements



a) Time Domain




b) Frequency Domain




# Time domain : Time Correlated Single Photon Counting



Electronic card to measure the delay between two pulses




- CFD : constant fraction discriminator
- TAC : Time to amplitude converter
- ADC : analog to digital converter






# Frequency domain : time / modulation detection



## Frequency domain : time / modulation detection



For a certain excitation frequency  $\omega$ , we can show


$$\tan \phi = \omega \tau_{\phi}$$

$\tau$  can be determined as well from the ratio of modulations  $m$  :

$$m = \frac{1}{\sqrt{1 + (\omega \tau_m)^2}}$$


# Time domain vs frequency domain

a) Time Domain



→ Useful for low light levels (single photon counting) : single molecules


b) Frequency Domain





→ Very fast measurements if only one frequency is used (average lifetime)

The fluorescence microscope

- Like in a spectrofluorimeter, it is necessary to separate excitation and emission wavelengths, as well as optical paths





$$\text{Resolution}_{x,y} = \lambda / 2[\eta \cdot \sin(\alpha)]$$


$$\text{Resolution}_z = 2\lambda / [\eta \cdot \sin(\alpha)]^2$$

Where  $\eta \cdot \sin(\alpha) = \text{NA}$  (Numerical aperture)


Example :  $\lambda = 550\text{nm}$ ,  $\text{NA} = 1.4$      $R_{xy} = 305\text{nm}$ ,  $R_z = 560\text{nm}$

## Illumination and collection modes


- EPISCOPIC ILLUMINATION
- CONFOCAL MICROSCOPY
- TOTAL INTERNAL REFLEXION MICROSCOPY
- MULTIPHOTON MICROSCOPY










## Wide field illumination



# ILLUMINATION EPISCOPIQUE





In a thick sample, light is collected from all sample planes




How to reject out-of-focus light ?

## Illumination and collection modes

- EPISCOPIC ILLUMINATION
- CONFOCAL MICROSCOPY
- TOTAL INTERNAL REFLEXION MICROSCOPY
- MULTIPHOTON MICROSCOPY




Then, to get an image, you need to do point scanning




Solution 1 : Laser is fixed, sample is moving

Then, to get an image, you need to do point scanning



**Epifluorescence**



**Section confocale**

